


Mr1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q95460

Item Number: CM0030714



Introduction

Rabbit polyclonal antibody against human MR1 (Q95460), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for MAIT cell and cancer immune surveillance research.

[Learn More](#)

Parameter	Value
Product Name	MR1 Rabbit pAb
Remarks/Alias	HLALS; MR1
Species	Human
GeneID (Human)	3140
GeneID	3140
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-291 of human MR1 (NP_001522.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q95460
Protein Weight	39kDa
Shipping	Ice bag