


Mrpl1 Rabbit Pab - Q9Byd6

Item Number: CM0024917



Introduction

Polyclonal antibody targeting human MRPL1 (L1mt, MRP-L1) for WB and ELISA. Detects the 39S ribosomal protein L1 in mitochondria, with cross-reactivity to human and rat. Suitable for mitochondrial translation studies.

[Learn More](#)

Parameter	Value
Product Name	MRPL1 Rabbit pAb
Remarks/Alias	L1MT; BM022; MRP-L1; MRPL1
Species	Human
GeneID (human)	65008
GeneID	65008
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 51-325 of human MRPL1 (NP_064621.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9BYD6
Protein Weight	37kDa
Shipping	Ice bag