


Mrpl20 Rabbit Monoclonal Antibody (Kd Validated) For Wb, IHC-P, Elisa - Q9Byc9

Item Number: CM0004205



Introduction

KD validated recombinant rabbit monoclonal antibody targeting mitochondrial ribosomal protein MRPL20. Suitable for Western blot, IHC-P, ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for mitochondrial translation studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MRPL20 Rabbit mAb
Remarks/Alias	L20mt; bL20m; MRP-L20
Species	Human
Gene ID	55052
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 48-147 of human MRPL20 (NP_060441.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYC9
Protein Weight	17kDa
Shipping	Ice bag