


Mrpl34 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bq48

Item Number: CM0023035



Introduction

Rabbit polyclonal antibody against MRPL34, a mitochondrial large ribosomal subunit protein. Validated for WB, IF/ICC, ELISA with cross-reactivity to human, mouse, and rat. Suitable for mitochondrial protein synthesis studies.

[Learn More](#)

Parameter	Specification
Product Name	MRPL34 Rabbit pAb
Remarks/Alias	L34mt; MRPL34
Species	Human
GeneID (Human)	64981
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-92 of human MRPL34 (NP_076426.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BQ48
Protein Weight	10kDa
Shipping	Ice bag