


Mrpl45 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Brj2

Item Number: CM0017283



Introduction

MRPL45 Rabbit Polyclonal Antibody, validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Targets mitochondrial large ribosomal subunit protein mL45 for ribosome research. Reliable IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	MRPL45 Rabbit pAb
Remarks/Alias	Mba1; L45mt; MRP-L45; MRPL45
Species	Human
GeneID (Human)	84311
GeneID	84311
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BRJ2
Protein Weight	35kDa
Shipping	Ice bag