


Mrps27 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92552

Item Number: CM0014284



Introduction

High-quality rabbit polyclonal antibody against human MRPS27 (mitochondrial ribosomal protein S27) validated for WB, IF/ICC, ELISA. Suitable for human, mouse, rat samples. Research use in mitochondrial translation and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	MRPS27 Rabbit pAb
Remarks/Alias	S27mt; MRP-S27; MRPS27
Species	Human
GeneID (Human)	23107
GeneID	23107
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 255-414 of human MRPS27 (NP_055899.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92552
Protein Weight	48kDa
Shipping	Ice bag