

Mrps36 Rabbit Pab For Wb, Ihc-P, Elisa - P82909

Item Number: CM0024671

Introduction

High-quality rabbit polyclonal antibody against human MRPS36/Alpha-ketoglutarate dehydrogenase component 4, validated for WB, IHC-P, and ELISA. Suitable for Human, Rat samples. Ideal for mitochondrial metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	MRPS36 Rabbit pAb
Remarks/Alias	MRPS36; DC47; MRP-S36; mitochondrial ribosomal protein S36
Species	Human
GeneID (Human)	6183
GeneID (Rat)	92259
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 1-103 of human MRPS36 (NP_150597.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P82909
Protein Weight	11kDa
Shipping	Ice bag