


Msr1 Rabbit Mab For Wb, Elisa - P21757

Item Number: CM0007038



Introduction

MSR1 Rabbit mAb is a monoclonal antibody for Western blot and ELISA, detecting human, mouse, and rat macrophage scavenger receptor 1 (MSR1/CD204) at ~50 kDa. Ideal for studying foam cell formation, atherosclerosis, and innate immune defense.

[Learn More](#)

Parameter	Value
Product Name	MSR1 Rabbit mAb
Remarks/Alias	SRA; SR-A; CD204; SR-AI; phSR1; phSR2; SCARA1; SR-AII; SR-AIII; MSR1
Species	Human
GeneID (Human)	4481
GeneID	4481
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MSR1 (P21757)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21757
Protein Weight	50kDa
Shipping	Ice bag