


Msr1/Sra/Sr-A Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P21757

Item Number: CM0024234



Introduction

Rabbit polyclonal antibody targeting human MSR1 (CD204/scavenger receptor A). Reacts with human, mouse, and rat. Validated for WB, IF-P, and ELISA. Immunogen: recombinant protein (aa 90-274 of human MSR1). Useful for atherosclerosis and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	MSR1 Rabbit pAb
Remarks/Alias	SRA; SR-A; CD204; SR-AI; phSR1; phSR2; SCARA1; SR-AII; SR-AIII; MSR1
Species	Human
GeneID (Human)	4481
GeneID	4481
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-274 of human MSR1 (NP_619729.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21757
Protein Weight	50kDa
Shipping	Ice bag