


Mst4 Rabbit Monoclonal Antibody For Wb And Elisa - Q9P289

Item Number: CM0003387



Introduction

Recombinant rabbit monoclonal antibody against MST4 (STK26), validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Suitable for studying cell growth regulation, apoptosis, and autophagy.

[Learn More](#)

Parameter	Value
Product Name	MST4 Rabbit mAb
Remarks/Alias	MASK; MST4
Species	Human
GeneID (Human)	51765
GeneID	51765
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 317-416 of human MST4 (Q9P289).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P289
Protein Weight	47kDa
Shipping	Ice bag