

Mt-Co3 Rabbit Polyclonal Antibody - Mt-Co3 Rabbit Pab

Item Number: CM0025859

Introduction

Rabbit anti-MT-CO3 polyclonal antibody for WB and ELISA applications, reacting with mouse and rat species. This antibody is raised against a synthetic peptide corresponding to amino acids 1-100 of mouse MT-CO3 (NP_904334.1). Suitable for mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	MT-CO3 Rabbit pAb
Remarks/Alias	COX3; MT-CO3
Species	Mouse
Gene ID (Human)	4514
Gene ID	17710
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse MT-CO3 (NP_904334.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
Protein Weight	29kDa
Shipping	Ice bag