


Mt-Nd1 Rabbit Polyclonal Antibody (Anti-Nadh Dehydrogenase Subunit 1) For Wb, Ihc-P, If/Icc, Elisa - P03886

Item Number: CM0023472



Introduction

Rabbit polyclonal antibody against human MT-ND1 (NADH dehydrogenase subunit 1). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Target molecular weight 36kDa. Ideal for mitochondrial complex I research.

[Learn More](#)

Parameter	Value
Product Name	MT-ND1 Rabbit pAb
Remarks/Alias	MTND1; ND1; MT-ND1
Species	Human
GeneID (Human)	4535
GeneID	4535
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MT-ND1 (YP_003024026.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P03886
Protein Weight	36kDa
Shipping	Ice bag