


Mt-Nd1 Rabbit Polyclonal Antibody For Wb And Elisa - P03888

Item Number: CM0025788



Introduction

High-quality MT-ND1 Rabbit Polyclonal Antibody validated for WB and ELISA. Recognizes human and mouse NADH-ubiquinone oxidoreductase chain 1, a key mitochondrial Complex I subunit. Ideal for mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	MT-ND1 Rabbit pAb
Remarks/Alias	TrnI; MT-ND1
Species	Mouse
Gene ID (Human)	4535
Gene ID	17716
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P03888
Protein Weight	36kDa
Shipping	Ice bag