

Mt-Nd3 Rabbit Polyclonal Antibody - Mt-Nd3

Item Number: CM0030955

Introduction

Rabbit polyclonal antibody targeting MT-ND3, raised against a synthetic peptide within amino acids 1-100 of mouse MT-ND3. Validated for WB and ELISA in human, mouse, and rat. Predicted molecular weight 13kDa. Ideal for mitochondrial complex I research.

[Learn More](#)

Parameter	Specification
Product Name	MT-ND3 Rabbit pAb
Remarks/Alias	ND3; MT-ND3
Species	Mouse
Gene ID (Human)	4537
Gene ID (Mouse)	17718
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse MT-ND3 (NP_904335.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Protein Weight	13kDa
Shipping	Ice bag