

Mt-Nd4 Rabbit Polyclonal Antibody - Mt-Nd4

Item Number: CM0030179

Introduction

This rabbit polyclonal antibody specifically targets MT-ND4 (ND4), a core subunit of mitochondrial complex I. Suitable for Western blot and ELISA applications, it recognizes human, mouse, and rat samples. Ideal for mitochondrial function and disease research.

[Learn More](#)

Parameter	Value
Product Name	MT-ND4 Rabbit pAb
Remarks/Alias	ND4; MT-ND4
Species	Mouse
Gene ID (Human)	4538
Gene ID (Mouse)	17719
Immunogen	Synthetic peptide corresponding to amino acids 360-459 of mouse MT-ND4 (NP_904337.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Protein Weight	51kDa
Shipping	Ice bag