


Mtap Rabbit Monoclonal Antibody For Wb, Elisa - Q13126

Item Number: CM0006213



Introduction

Rabbit monoclonal antibody targeting human MTAP for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Suitable for studying methionine salvage and polyamine biosynthesis pathways.

[Learn More](#)

Parameter	Value
Product Name	MTAP Rabbit mAb
Remarks/Alias	BDMF; MSAP; DMSFH; LGMBF; DMSMFH; c86fus; HEL-249; MTAP
Species	Human
GeneID (Human)	4507
GeneID	4507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human MTAP (NP_002442.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13126
Protein Weight	31kDa
Shipping	Ice bag