

Mtf2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y483

Item Number: CM0025096



Introduction

Rabbit polyclonal antibody against human MTF2, validated for WB and ELISA. Reacts with human, mouse, rat. Suitable for PRC2 complex studies and ES cell differentiation research.

[Learn More](#)

Parameter	Value
Product Name	MTF2 Rabbit pAb
Remarks/Alias	M96; PCL2; TDRD19A; dj976O13.2; MTF2
Species	Human
GeneID (human)	22823
GeneID	22823
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human MTF2 (NP_001157864.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y483
Protein Weight	67kDa
Shipping	Ice bag