

Mthfr Rabbit Polyclonal Antibody For Wb, Elisa - P42898

Item Number: CM0018491



Introduction

Rabbit polyclonal antibody targeting human methylenetetrahydrofolate reductase (MTHFR), validated for Western blot and ELISA, with cross-reactivity to rat. Essential for studying folate and methionine metabolism, and related disease variants.

[Learn More](#)

Parameter	Value
Product Name	MTHFR Rabbit pAb
Remarks/Alias	MTHFR
Species	Human
GeneID (Human)	4524
GeneID	4524
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human MTHFR (NP_005948.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P42898
Protein Weight	75kDa
Shipping	Ice bag