


Mtif2 Rabbit Pab - P46199

Item Number: CM0023498



Introduction

Polyclonal antibody targeting mitochondrial translation initiation factor IF-2 (MTIF2), suitable for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studying mitochondrial protein synthesis and related pathways.

[Learn More](#)

Parameter	Value
Product Name	MTIF2 Rabbit pAb
Remarks/Alias	MTIF2
Species	Human
GeneID (Human)	4528
GeneID	4528
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 428-727 of human MTIF2 (NP_002444.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46199
Protein Weight	81kDa
Shipping	Ice bag