

## Mtif2 Rabbit Pab - P46199

Item Number: CM0023468

### Introduction

Polyclonal antibody targeting mitochondrial translation initiation factor IF-2 (MTIF2), suitable for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studying mitochondrial protein synthesis and related pathways.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | MTIF2 Rabbit pAb                                                                                                   |
| Remarks/Alias    | MTIF2                                                                                                              |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 4528                                                                                                               |
| GeneID           | 4528                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 428-727 of human MTIF2 (NP_002444.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | P46199                                                                                                             |
| Protein Weight   | 81kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |