

Mtor Rabbit Pab - P42345

Item Number: CM0019663



Introduction

mTOR Rabbit Polyclonal Antibody validated for WB, IF/ICC, IP, ELISA with human, mouse, rat reactivity. Central regulator of cell growth and metabolism; suitable for signaling and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	mTOR Rabbit pAb
Remarks/Alias	SKS; FRAP; FRAP1; FRAP2; RAFT1; RAPT1; mTOR
Species	Human
GeneID (Human)	2475
GeneID	2475
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human mTOR (NP_004949.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42345
Protein Weight	289kDa
Shipping	Ice bag