

Mtx1 Rabbit Polyclonal Antibody (Pab) For Wb, Elisa - Q13505

Item Number: CM0024701

Introduction

Rabbit polyclonal antibody targeting MTX1 (Metaxin-1), a mitochondrial outer membrane protein involved in protein import. Validated for WB and ELISA, cross-reactive with human, mouse, and rat. Ideal for mitochondrial transport research.

[Learn More](#)

Parameter	Value
Product Name	MTX1 Rabbit pAb
Remarks/Alias	MTX; MTXN; MTX1
Species	Human
GeneID (Human)	4580
GeneID	4580
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 186-435 of human MTX1 (NP_942584.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13505
Protein Weight	51kDa
Shipping	Ice bag