

Mx1 Rabbit Pab - P20591

Item Number: CM0021127



Introduction

Interferon-induced GTP-binding protein Mx1 polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human and mouse MX1. Ideal for antiviral innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	MX1 Rabbit pAb
Remarks/Alias	MX; MxA; IFI78; IFI-78K; IncMX1-215; MX1
Species	Human
GeneID (human)	4599
GeneID	4599
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of human MX1 (NP_001171517.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P20591
Protein Weight	76kDa
Shipping	Ice bag