

Myb Rabbit Pab For Wb, Elisa - P10242

Item Number: CM0017661

Introduction

MYB rabbit polyclonal antibody (P10242) for WB and ELISA. Derived from synthetic peptide immunogen within amino acids 380-640 of human MYB. Suitable for human samples in transcriptional research.

[Learn More](#)

Parameter	Value
Product Name	MYB Rabbit pAb
Remarks/Alias	efg, Cmyb, c-myb, c-myb_CDS, MYB
Species	Human
GeneID (Human)	4602
GeneID	4602
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 380-640 of human MYB (NP_005366.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P10242
Protein Weight	72kDa
Shipping	Ice bag