

Myf5 Rabbit Pab For Wb - P13349

Item Number: CM0022330



Introduction

MYF5 Rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and mouse MYF5, a transcription factor critical for muscle differentiation. Ideal for myogenesis research and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	MYF5 Rabbit pAb
Remarks/Alias	EORVA; bHLHc2; MYF5
Species	Human
Gene ID (Human)	4617
Gene ID	4617
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 143-220 of human MYF5 (NP_005584.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13349
Protein Weight	28kDa
Shipping	Ice bag