

Myf6 Rabbit Pab For Wb And Elisa - P23409

Item Number: CM0013270

Introduction

The MYF6 Rabbit pAb (P23409) is a polyclonal antibody for Western blot and ELISA applications. It detects human, mouse, and rat MYF6, a 27 kDa transcription factor critical for skeletal muscle differentiation.

[Learn More](#)

Parameter	Value
Product Name	MYF6 Rabbit pAb
Remarks/Alias	CNM3; MRF4; myf-6; bHLHc4; MYF6
Species	Human
GeneID (human)	4618
GeneID	4618
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-180 of human MYF6 (NP_002460.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23409
Protein Weight	27kDa
Shipping	Ice bag