


Myh14 Rabbit Mab For Wb, If-P, Elisa - Q7Z406

Item Number: CM0004360



Introduction

MYH14 Rabbit monoclonal antibody validated for WB, IF-P, ELISA; reactive with human, mouse, rat. Suitable for studying cytokinesis, cell shape, and non-muscle myosin functions.

[Learn More](#)

Parameter	Value
Product Name	MYH14 Rabbit mAb
Remarks/Alias	DFNA4; MHC16; MYH17; PNMHH; DFNA4A; myosin; FP17425; NMHC II-C; NMHC-II-C; MYH14
Species	Human
GeneID (Human)	79784
GeneID	79784
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 600-700 of human MYH14 (NP_079005.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z406
Protein Weight	228kDa
Shipping	Ice bag