

Myh6 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P13533

Item Number: CM0010574

Introduction

Rabbit monoclonal antibody against human MYH6 (alpha-MHC) suitable for WB, IF-P, ELISA; cross-reacts with mouse. For cardiac muscle contraction and cardiomyopathy-associated mutation studies. SWISS: P13533.

[Learn More](#)

Parameter	Value
Product Name	MYH6 Rabbit mAb
Remarks/Alias	ASD3; MYHC; SSS3; CMH14; MYHCA; CMD1EE; alpha-MHC
Species	Human
GeneID (Human)	4624
GeneID	4624
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 160-260 of human MYH6 (NP_002462.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13533
Protein Weight	224kDa
Shipping	Ice bag