

Myh6/A-Mhc Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P13533

Item Number: CM0018022

Introduction

High-quality rabbit polyclonal antibody against MYH6/α-MHC, validated for WB, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for cardiac muscle research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	MYH6/α-MHC Rabbit pAb
Remarks/Alias	ASD3, MYHC, SSS3, CMH14, MYHCA, CMD1EE, alpha-MHC, MYH6/α-MHC
Species	Human
Gene ID (Human)	4624
Gene ID	4624
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 600-700 of human MYH6/α-MHC (NP_002462.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P13533
Protein Weight	224 kDa
Shipping	Ice bag