

Myh6/A α -Mhc Rabbit Mab - Q02566

Item Number: CM0010348



Introduction

Anti-MYH6/A α -MHC rabbit monoclonal antibody, validated for WB, IF-P, IHC-P, ELISA. Detects mouse and rat myosin-6 (alpha cardiac muscle myosin heavy chain). Ideal for cardiac contractile research and muscle biology studies. SWISS: Q02566.

[Learn More](#)

Parameter	Value
Product Name	MYH6/A α -MHC Rabbit mAb
Remarks/Alias	Myhca; Myhc-a; alphaMHC; alpha-MHC; A830009F23Rik
Species	Mouse
GeneID (human)	4624
GeneID	17888
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 621-720 of mouse MYH6/A α -MHC (NP_034986.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q02566
Protein Weight	224 kDa
Shipping	Ice bag