

Myh6/Myh7 Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - P13533 / P12883

Item Number: CM0009406

Introduction

Recombinant rabbit PolymAb® against human MYH6/MYH7 (cardiac myosin heavy chain alpha/beta). Validated for WB, IF-P, IHC-P, ELISA with human, mouse, rat reactivity. Ideal for cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	MYH6/MYH7 Rabbit PolymAb®
Remarks/Alias	ASD3; MYHC; SSS3; CMH14; MYHCA; CMD1EE; alpha-MHC;CMH1; MPD1; SPMD; SPMM; CMD1S; MYHCB
Species	Human
GeneID (human)	4624/4625
GeneID	4624/4625
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1730-1830 of human MYH6/MYH7 (NP_002462.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13533/P12883
Protein Weight	223kDa
Shipping	Ice bag