

MyL2 Rabbit Mab For Wb, If-P, Elisa - P10916

Item Number: CM0011507



Introduction

MYL2 Rabbit mAb is a monoclonal antibody validated for WB, IF-P, ELISA. It targets human MYL2 (MLC2) and cross-reacts with mouse and rat. Ideal for cardiac muscle research.

[Learn More](#)

Parameter	Value
Product Name	MYL2 Rabbit mAb
Remarks/Alias	MLC2; CMH10; MFM12; MLC-2; MLC-2v; MLC-2s/v; MYL2
Species	Human
Gene ID (Human)	4633
Gene ID	4633
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P10916
Protein Weight	19kDa
Shipping	Ice bag