

MyL6B Rabbit Pab For Wb, Elisa - P14649

Item Number: CM0019565



Introduction

Rabbit polyclonal antibody against human MYL6B (myosin light chain 6B), validated for WB and ELISA; cross-reacts with mouse and rat; suitable for muscle biology and myosin research.

[Learn More](#)

Parameter	Value
Product Name	MYL6B Rabbit pAb
Remarks/Alias	MLC15A; MYL6B
Species	Human
GeneID (Human)	140465
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-208 of human MYL6B (NP_002466.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14649
Protein Weight	23kDa
Shipping	Ice bag