

Mylk Rabbit Mab For Wb, Elisa - Q15746

Item Number: CM0008691



Introduction

Rabbit monoclonal antibody against human MYLK (MLCK) for WB and ELISA. Suitable for smooth muscle contraction, cell migration, and inflammatory response research. Molecular weight 211 kDa.

[Learn More](#)

Parameter	Value
Product Name	MYLK Rabbit mAb
Remarks/Alias	KRP; AAT7; MLCK; MLCK1; MMIHS; MYLK1; MMIHS1; MYLK-L; smMLCK; MLCK108; MLCK210; MSTP083; MYLK
Species	Human
GeneID (Human)	4638
GeneID	4638
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1800-1900 of human MYLK (Q15746).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15746
Protein Weight	211kDa
Shipping	Ice bag