

Mylk3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q32Mk0

Item Number: CM0026703

Introduction

Rabbit polyclonal antibody targeting MYLK3 (Cardiac-MLCK), suitable for WB, IF/ICC, ELISA; cross-reacts with human and mouse, detects ~88 kDa protein. Ideal for cardiac sarcomere and kinase signaling studies.

[Learn More](#)

Parameter	Value
Product Name	MYLK3 Rabbit pAb
Remarks/Alias	MLCK; MLCK2; caMLCK; MYLK3
Species	Human
GeneID (human)	91807
GeneID	91807
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human MYLK3 (NP_872299.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q32MK0
Protein Weight	88kDa
Shipping	Ice bag