

Myo10 Rabbit Pab For Wb, Elisa - Q9Hd67

Item Number: CM0018258



Introduction

MYO10 Rabbit pAb (Q9HD67) is a polyclonal antibody validated for WB and ELISA, targeting human MYO10 with cross-reactivity to mouse. Ideal for studying myosin-X functions in cell adhesion, filopodia, and intracellular transport.

[Learn More](#)

Parameter	Value
Product Name	MYO10 Rabbit pAb
Remarks/Alias	MyoX; MYO10
Species	Human
Gene ID	4651 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 845-944 of human MYO10 (NP_036466.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HD67
Protein Weight	237kDa
Shipping	Ice bag