

Myo5A Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y4I1

Item Number: CM0025685



Introduction

MYO5A Rabbit pAb is a polyclonal antibody raised in rabbit against recombinant human MYO5A protein. Validated for WB and ELISA in human, mouse, and rat. Suitable for research on myosin-Va motor function, melanosome transport, and neuronal polarization.

[Learn More](#)

Parameter	Value
Product Name	MYO5A Rabbit pAb
Remarks/Alias	GS1; MYO5; MYH12; MYR12; MYO5A
Species	Human
GeneID (Human)	4644
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 950-1150 of human MYO5A (NP_000250.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4I1
Protein Weight	215kDa
Shipping	Ice bag