


Myoc Rabbit Pab For Wb, Elisa - Q99972

Item Number: CM0024287



Introduction

High-quality rabbit polyclonal antibody targeting human MYOC (Myocilin). Validated for WB and ELISA applications, with cross-reactivity to human and mouse. Ideal for glaucoma, cell adhesion, and migration research.

[Learn More](#)

Parameter	Value
Product Name	MYOC Rabbit pAb
Remarks/Alias	GPOA; JOAG; TIGR; GLC1A; JOAG1; MYOC
Species	Human
GeneID (Human)	4653
GeneID	4653
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99972
Protein Weight	57kDa
Shipping	Ice bag