

Myod1 Rabbit Pab For Wb And Elisa - P15172

Item Number: CM0016729

Introduction

MyoD1 Rabbit Polyclonal Antibody targeting human MyoD1 (SWISS: P15172) provides reliable detection in WB and ELISA with cross-reactivity to mouse, supporting muscle differentiation and myogenesis research.

[Learn More](#)

Parameter	Value
Product Name	MyoD1 Rabbit pAb
Remarks/Alias	PUM; MYF3; MYOD; CMYP17; bHLHc1; MYODRIF; MyoD1
Species	Human
Gene ID (Human)	4654
Gene ID	4654
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 226-300 of human MyoD1 (NP_002469.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P15172
Protein Weight	35kDa
Shipping	Ice bag