


Myoglobin Rabbit Monoclonal Antibody For Wb And Elisa - P02144

Item Number: CM0003517



Introduction

Rabbit monoclonal antibody targeting human Myoglobin (P02144), suitable for WB and ELISA applications. Detects mouse and rat Myoglobin. Essential for studies on oxygen transport, muscle physiology, and nitrite reductase activity.

[Learn More](#)

Parameter	Value
Product Name	Myoglobin Rabbit mAb
Remarks/Alias	MYOSB; PVALB; Myoglobin
Species	Human
GeneID (Human)	4151
GeneID	4151
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 55-154 of human Myoglobin (P02144).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P02144
Protein Weight	17kDa
Shipping	Ice bag