

## Myom2 Rabbit Polyclonal Antibody For Wb, Elisa - P54296

Item Number: CM0019029

### Introduction

MYOM2 Rabbit pAb is a polyclonal antibody raised against human MYOM2. Validated for WB and ELISA, it cross-reacts with mouse and rat, suitable for studying myofibrillar M band architecture and sarcomere organization.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | MYOM2 Rabbit pAb                                                                                                 |
| Remarks/Alias    | TTNAP; MYOM2                                                                                                     |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 9172                                                                                                             |
| GeneID           | 9172                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human MYOM2 (NP_003961.3) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | P54296                                                                                                           |
| Protein Weight   | 165kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                          |