


Myot Rabbit Pab For Wb, If/Icc, Elisa - Q9Ubf9

Item Number: CM0026831



Introduction

High-quality rabbit polyclonal antibody targeting human myotilin (MYOT). Validated for WB, IF/ICC, and ELISA applications with cross-reactivity to mouse and rat. Ideal for muscle cell and myofibril assembly research.

[Learn More](#)

Parameter	Value
Product Name	MYOT Rabbit pAb
Remarks/Alias	MFM3; TTID; TTOD; LGMD1; LGMD1A; MYOT
Species	Human
GeneID (Human)	9499
GeneID	9499
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 259-314 of human MYOT (NP_001129412.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBF9
Protein Weight	55kDa
Shipping	Ice bag