

Myoz2 Rabbit Pab For Wb, If/Icc, Elisa - Q9Npc6

Item Number: CM0028597

Introduction

MYOZ2 Rabbit pAb is a polyclonal antibody applicable in WB, IF/ICC, ELISA, recognizing human, mouse, rat Myozenin-2 (calsarcin-1). Key in sarcomere Z-line organization and calcineurin signaling, suitable for cardiac/skeletal muscle studies.

[Learn More](#)

Parameter	Value
Product Name	MYOZ2 Rabbit pAb
Remarks/Alias	CS-1; CMH16; C4orf5; FATZ-2; MYOZ2
Species	Human
GeneID (Human)	51778
GeneID	51778
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPC6
Protein Weight	30kDa
Shipping	Ice bag