

N5-Methylcytidine Rabbit Pab - M5C

Item Number: CM0024615

Introduction

Rabbit polyclonal antibody targeting N5-Methylcytidine (m5C), a key RNA modification. Validated for IF/ICC and ELISA with human cross-reactivity, supporting epitranscriptomic research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	N5-Methylcytidine Rabbit pAb
Remarks/Alias	m5C; N5-Methylcytidine
Immunogen	Chemical compounds corresponding to N5-Methylcytidine / m5C.
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
Shipping	Ice bag