

N6-Methyladenosine / M6A Mouse Mab - A19136

Item Number: CM0003140

Introduction

N6-methyladenosine (m6A) mouse monoclonal antibody validated for dot blot (DB) and ELISA. Species-independent detection of m6A RNA modification. Ideal for epigenetic and RNA research applications.

[Learn More](#)

Parameter	Value
Product Name	N6-methyladenosine / m6A Mouse mAb
Remarks/Alias	N6-methyladenosine; m6A; N6-methyladenosine / m6A
Gene ID (Human)	m6A
Immunogen	Chemical compounds corresponding to N6-methyladenosine / m6A
Source	Mouse
Category	Monoclonal Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag