

## N6-Methyladenosine (M6A) Rabbit Monoclonal Antibody For Merip, Elisa, If/Icc, Db And Nucleotide Array

Item Number: CM0000064

### Introduction

High-quality N6-methyladenosine (m6A) Rabbit mAb for meRIP, ELISA, IF/ICC, DB, and nucleotide array applications. Exhibits species-independent cross reactivity. Suitable for epitranscriptomic studies and RNA modification mapping.

[Learn More](#)

| Parameter        | Value                                                        |
|------------------|--------------------------------------------------------------|
| Product Name     | N6-methyladenosine / m6A Rabbit mAb                          |
| Remarks/Alias    | N6-methyladenosine; m6A; N6-methyladenosine / m6A            |
| Gene ID (Human)  | m6A                                                          |
| Immunogen        | Chemical compounds corresponding to N6-methyladenosine / m6A |
| Source           | Rabbit                                                       |
| Category         | Other Modified Antibodies                                    |
| Application      | DB, IF/ICC, ELISA, meRIP, Nucleotide Array                   |
| Cross Reactivity | Species independent                                          |
| Shipping         | Ice bag                                                      |