

N6-Methyladenosine (M6A) Rabbit Polyclonal Antibody For Merip, Dot Blot, Elisa, And Nucleotide Array

Item Number: CM0000050

Introduction

Species-independent anti-m6A rabbit polyclonal antibody validated for meRIP, dot blot, ELISA, IF/ICC, and nucleotide array. Suitable for RNA modification research and epitranscriptomic assay development.

[Learn More](#)

Parameter	Value
Product Name	N6-methyladenosine / m6A Rabbit pAb
Remarks / Alias	N6-methyladenosine; m6A; N6-methyladenosine / m6A
Gene ID (Human)	m6A
Immunogen	Chemical compounds corresponding to N6-methyladenosine / m6A
Source	Rabbit
Category	Other Modified Antibodies
Application	DB, IF/ICC, ELISA, meRIP, Nucleotide Array
Cross Reactivity	Species independent
Shipping	Ice bag