


Naa60 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H7X0

Item Number: CM0028487



Introduction

NAA60 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Reacts with Human, Mouse, Rat. N-alpha-acetyltransferase 60 (NatF) acetylates transmembrane proteins, essential for chromosome segregation. Ideal for cancer and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	NAA60 Rabbit pAb
Remarks/Alias	HAT4; NatF; NAT15; hNaa60; NAA60
Species	Human
GeneID (Human)	79903
GeneID	79903
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-184 of human NAA60 (NP_001077070.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H7X0
Protein Weight	27kDa
Shipping	Ice bag