


Naf1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Hr8

Item Number: CM0022331



Introduction

Generated against recombinant fusion protein corresponding to amino acids 150-400 of human NAF1, this antibody recognizes human, mouse, rat NAF1 for Western blot and ELISA applications in ribosome biogenesis research.

[Learn More](#)

Parameter	Value
Product Name	NAF1 Rabbit pAb
Remarks/Alias	NAF1
Species	Human
GeneID	92345
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-400 of human NAF1 (NP_612395.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96HR8
Protein Weight	54kDa
Shipping	Ice bag