


Nagk Rabbit Pab For Wb, Elisa - Q9Uj70

Item Number: CM0030282



Introduction

NAGK rabbit polyclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studies on GlcNAc kinase function in glycan metabolism and innate immunity.

[Learn More](#)

Parameter	Specification
Product Name	NAGK Rabbit pAb
Remarks/Alias	GNK; HSA242910; NAGK
Species	Human
GeneID (Human)	55577
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-210 of human NAGK (NP_060037.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJ70
Protein Weight	37kDa
Shipping	Ice bag